



Washers & Sanitizers to automatically wash and apply disinfectants to wheelchairs and other hard-to-clean medical equipment.

AQ-8200

Features & Options

DESIGNED TO ACCOMODATE MOST ENDOSCOPE CASSETTES

Standard Features:

- Double Wall Stainless Steel Chamber
- Stainless Steel Spray Arms
- Carousel Style Loading & Wash System
- Adjustable Time Cycles
- Wash and Rinse Cycles
- Automatic Chemical Injection
- LCD Status Display
- Front Door Window
- Adjustable Leveling Feet
- Fill & Drain Hose Set
- 115V, 50/60Hz, 20 Amp AC Operation



Optional Features:

- Pass Through Option-2 Doors (Shown)
- Door Interlock System
- Disinfectant Application System
- 230V, 50/60Hz, Single Phase AC Operation



Total Outer Dimensions

- **Width** 36"
- **Length** 31-5/8"
- **Height** 58-1/2"

Effective Chamber Dimensions

- **Width** 27-1/4"
- **Length** 27-1/4"
- **Height** 26-1/2"

Like all Aqua Phase models, in order to operate the machine once the unit is located, all you need to do is hook up two hoses — one for the hot water fill and one for the drain, connect the Aqua Phase to an appropriate power source and begin saving time, improving efficiency and enhancing your infection control protocol.